

INSTALLATION GUIDE

1. Create the folder C:\MES
2. Download MES6.5_released.zip from the above link and unzip it to the newly created MES folder
3. Start MATLAB as an administrator
4. Use the MATLAB method 'Set Path', add the folder C:\MES\MES6.5_released\main, and save this new path
5. Type 'mes' into the MATLAB command window
6. Type Y to create the 'Local Storage' directory
7. Two default MES user accounts will be created: admin and service
8. Use 'admin' for work. By default its password is empty and will not be required to input
9. From the second startup, you no more need to run MATLAB as an administrator
10. The demo license in the package contains the 'main', 'feldolg' and 'analysis' modules, and is suitable for viewing the sample measurement file and provides specialized line-scan, image and curve analysis software elements.
If you already have a license containing more modules, you may use it with this package by copying it to the 'settings' folder.
To request access for other modules, please use our License request application.